

APC 3

Argon plasma coagulation module for the VIO® 3



APC 3

Argon plasma coagulation module for the VIO® 3

- Operation via VIO® 3
- Plug connection of VIO® 3 with integrated Erbe Communication Bus (ECB)
- 1 APC socket exclusively for FiAPC® instruments from Erbe
- Extended VIO® 3: 5 HF instruments plus 1 APC instrument
- Gas supply via gas cylinder or central building supply
- Argon gas with 4.8 level of purity
- Min. gas pressure 5 bar
- Modes:
 - forcedAPC
 - pulsedAPC
 - preciseAPC



APC socket

Filter integrated APC connector (FiAPC®):

- Patented
- Filters particles from the gas cylinder
- Prevents return flow

Instrument recognition
(plug-and-play):

- Mode
- Effect
- Flow



APC 3 front

Plug connection to VIO® 3:
power and ECB communication

