

No: 10340-000

ESM 2 Suction module



- The type of activation and the degree of suction of the ESM 2 are adjusted via the interface of the ERBEJET 2
- Activation occurs either simultaneous to the application of the waterjet
- The suction pump can be programmed to continue running for a certain time after waterjet is no longer activated.
- Good visibility of the target tissue

Mains connection

System voltage	120 - 240 Volt
Power frequency	50 / 60 Hz

Additional product information

Communication Bus	{ERBE Communication Bus (ECB)}
Packaging unit	1 Piece

Standards

Protection category acc. to EN 60 601-1	I
Classification acc. to 93/42 EWG	IIa
Type acc. to EN 60 601-1	CF