

No: 10402-000

ERBECRYO[®] 2



- Convenient operation with plug and play: ERBECRYO 2 automatically sets the appropriate parameters for every probe
- Even more convenience: all important settings appear on the display: - Effect levels - Timer function
- Activation using the footswitch
- Reproducible freeze performance from the first time the probe is used through to the last

Mains connection

| | |
|-----------------------------------|------------------------|
| System voltage | 100-120 / 220-240 Volt |
| Power frequency | 50 / 60 Hz |
| Potential equalization connection | yes |

Dimensions and weight

| | |
|------------------------|--------------------|
| Width x Height x Depth | 410 X 130 X 370 mm |
| Weight | 6,7 kg |

Additional product information

| | |
|------------------------|---------|
| Packaging unit | 1 Piece |
| Programmspeicherplätze | {10} |

Standards

| | |
|---|-----|
| Classification acc. to 93/42 EWG | IIb |
| Protection category acc. to EN 60 601-1 | I |
| Type acc. to EN 60 601-1 | CF |